

# Manual & Warranty Card


# ART TRAY

Bite tray for real impression

# COVER

This 'Tray Cover' is used by putting it on the 'ART Tray' when taking dental impressions.

## ■ Product Description



## ■ Product Features

- Two layers of nonwoven are attached to one
- Harmless material

## ■ Product Precautions

- Store in a cool place away from moisture and heat
- Use proper force to disconnect the product from the connected area
- When inserting the 'ART tray' into the cover, make sure that the frame direction and cover direction are correct

## ■ Product Warranty

- This product is guaranteed for one year from the date of purchase.
- After purchase, if abnormality is found or problems occurred, prepare this warranty card and transaction receipt, and contact the place of purchase to request for repair or replacement.

## ■ How to use

Prepare the 'cover paper'	Carefully hold the connected area and remove the 'cover'	When the separation is completed, gently rubs the insert part to make a gap	Check the correct direction and insert it into the 'ART Tray'
---------------------------	----------------------------------------------------------	-----------------------------------------------------------------------------	---------------------------------------------------------------



ART dental lab  
FIXED - REMOVABLE - IMPLANTS

703-378-8555  
artdlab1@gmail.com


# Manual & Warranty Card

# ART TRAY

Bite tray for real impression

'ART Tray' is the bite tray for taking dental impressions

## ■ Product Description



## ■ Product Features

- Reusable & Eco-friendly
- Non-toxic stainless steel material
- More accurate impression work
- Easy to remove residues after work
- Available in high-temperature, high-pressure sterilizers (Autoclave)


## ■ Product Precautions

- Use 'ART Tray' only for original purpose
- Sterilize and disinfect before first use
- When inserting the cover, bend the wings with proper force

## ■ Product Warranty

- This product is guaranteed for one year from the date of purchase.
- After purchase, if abnormality is found or problems occurred, prepare this warranty card and transaction receipt, and contact the place of purchase to request for repair or replacement.

## ■ How to use



Hold the handle and get ready

Bend wings slightly

Insert the cover in the proper direction

Push it in completely and you're ready